**Gewissheit in 15 min**

# **AESKU.**RAPID

SARS-CoV-2 Rapid Test

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# Evaluating the clinical utility and sensitivity of SARS-CoV-2 antigen testing in relation to RT-PCR Ct values

Lukas Lanser<sup>1</sup> · Rosa Bellmann-Weiler<sup>1</sup> · Karla-Wanda Öttl<sup>1</sup> · Lukas Huber<sup>1</sup> · Andrea Griesmacher<sup>2</sup> · Igor Theurl<sup>1</sup> · Günter Weiss<sup>1</sup>

Received: 4 October 2020 / Accepted: 21 October 2020  
© Springer-Verlag GmbH Germany, part of Springer Nature 2020

To the Editor,

Diagnosis of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) in patients with suspected coronavirus disease 2019 (COVID-19) is most widely performed with real time polymerase chain reaction (RT-PCR) considered as gold standard [1]. RT-PCR is a highly sensitive diagnostic method for detecting viral ribonucleic acid (RNA) with the disadvantages of logistics for transport of samples to specific laboratories and the relative long duration of the diagnostic method. Positive SARS-CoV-2 RT-PCR also do not allow definitive conclusions whether the subject is still contagious or not. This can be partly attained by establishing the cycle threshold (Ct) value depicting the particular amount of viral RNA in the sample and thus allowing a conclusion on the viral load and infectivity [2]. However, because of the diagnostic effort and duration of the test, faster and less laborious tests attract interest which may help to rapidly identify and contain infected individuals. Therefore, point-of-care-testing with antigen tests, providing results in a couple of minutes, is currently evaluated for routine clinical use.

We compared the SARS-CoV-2 antigen detection in nasopharyngeal swab samples by the Panbio™ COVID-19 Ag Rapid test (Abbott, Chicago, Illionis) with the simultaneous routinely conducted RT-PCR analysis of SARS-CoV-2 orf1 RNA detection with the cobas® analyzer (Roche Diagnostics GmbH, Mannheim, Germany). The nasopharyngeal swab samples were collected from 53 patients with PCR-confirmed SARS-CoV-2 infection during their hospital stay in different stages of the disease. Panbio™ COVID-19 Ag

Rapid test was performed right after nasopharyngeal swab sampling while RT-PCR was routinely performed in our central laboratory facility. RT-PCR was negative in two patients suggesting an already subsided infection; consistent with it the Panbio™ COVID-19 Ag Rapid test was also negative. Among 51 RT-PCR SARS-CoV-2 positive patients, the Panbio™ COVID-19 Ag Rapid test was positive in 31 subjects depicting a poor sensitivity of 60.8% (95% CI 46.1–74.2%), compared to 93.3% in the manufacturer's information. In the 14 patients with a Ct-value  $\leq 25$ , being indicative for higher viral loads, the sensitivity for the Panbio™ COVID-19 Ag Rapid test was at a level of 85.7% (95% CI 57.2–98.2%, Table 1). Panbio™ COVID-19 Ag Rapid test was positive in 36.4% (95% CI 17.2–59.3%) of the patients with a Ct-value  $> 30$  (Table 1). Of note, when we included subjects with a Ct value  $\leq 30$ , which is considered to be a threshold for infectivity and currently recommended by the German Robert Koch Institute as an important cutoff to identify SARS-CoV2 contagious subjects, we found that the Ag rapid test correctly identified 79.3% of individuals. However, looking on the other side of the coin we found that subjects with a Ct-value  $> 30$  were antigen positive in 36.4% of cases. The test is still positive in a considerable number of patients which are considered as being non-infectious according to the German Robert Koch Institute. Changing the cut-off value to Ct  $> 33$  reduced the positive results of the Panbio™ COVID-19 Ag Rapid test to 16.7% (Table 1) whereas no positive Panbio™ COVID-19 Ag Rapid test result was achieved in subjects with RT-PCR Ct-values  $\geq 35$  ( $n = 4$ ).

These results indicate a poor performance of the Panbio™ COVID-19 Ag Rapid test detecting SARS-CoV-2 antigen in patients with low virus load but a good performance in patients with high virus load. Viral load of SARS-CoV-2 was reported to be the highest around the time of symptom onset, and most probably become undetectable within approximately 2 weeks. Additionally, patients with more severe symptoms also seems to have a higher viral

✉ Günter Weiss  
guenter.weiss@i-med.ac.at

<sup>1</sup> Department of Internal Medicine II, Innsbruck Medical University, Anichstraße 35, 6020 Innsbruck, Austria

<sup>2</sup> Central Institute for Medical and Chemical Laboratory Diagnosis, Innsbruck University Hospital, Innsbruck, Austria

**Table 1** Antigen test results stratified by different Ct value cut-offs

	Antigen test negative <i>n</i> (% within Ct value)	Antigen test positive <i>n</i> (% within Ct value)	Total <i>n</i> (% within total)
Ct value > 25	18 (48.6%)	19 [51.4% (95% CI 34.4–68.1)]	37 (72.5%)
Ct value ≤ 25	2 (14.3%)	12 [85.7% (95% CI 57.2–98.2)]	14 (27.5%)
Ct value > 30	14 (63.6%)	8 [36.4% (95% CI 17.2–59.3)]	22 (43.1%)
Ct value ≤ 30	6 (20.7%)	23 [79.3% (95% CI 60.3–92.0)]	29 (56.9%)
Ct value > 33	5 (83.3%)	1 [16.7% (95% CI 0.42–64.1)]	6 (11.8%)
Ct value ≤ 33	15 (33.3%)	30 [66.7% (95% CI 51.1–80.0)]	45 (88.2%)
Total	20 (39.2%)	31 [60.8% (95% CI 46.1–74.2)]	51 (100%)

loads [3, 4]. Therefore, antigen tests, like the one investigated herein, might be suitable and effective for rapidly identifying infectious subjects with symptoms compatible to a COVID-19 infection in a primary care setting given its sensitivity of around 80%. This is of relevance specifically in epidemic situations with numerous other circulating virus such as influenza but also for rapidly identifying SARS-CoV-2 positive patients for immediate initiation of containment and contact tracing strategies for better controlling the spread of the infection mainly during the cold season [5]. Moreover, the test may be of specific benefit in subjects with short onset of symptoms as viral loads are highest during the first days of infection and because the test has been shown herein to have an even higher sensitivity in patients with higher viral loads. However, one has to admit that not all subjects with a Ct-value ≤ 30 have been correctly identified as being SARS-CoV-2 infected. This would warrant the development of additional strategies to minimize the risk of false negative results. Such algorithms could include the performance of a RT-PCR test in subjects with symptoms compatible with COVID-19 but a negative antigen result. The discrepancy between the real sensitivity and that claimed by the manufacturer of SARS-CoV-2 Ag tests was already demonstrated recently [6]. Since these are CE-marked diagnostic devices this question the reliability of the European admission process.

The low number of positive antigen test results in subjects with higher Ct value > 33 which are considered to be not contagious (<https://www.cdc.gov/coronavirus/2019-ncov/hcp/duration-isolation.html>) could have implications for follow up screening of initially antigen positive subjects to determine the time point when they are no longer contagious, which is specifically relevant for estimating the duration of quarantine measures. Whether or not this antigen test is suitable for public screening measures mostly of asymptomatic patients to identify undetected subjects with COVID-19 without generating a considerable number of false positive results remains to be shown. Although our results are quite encouraging regarding the use of antigen test as point-of-care diagnostic which may contribute to a better control of the pandemic, we need

longitudinal studied which larger patient cohorts to corroborate our results and to develop algorithms or to identify those subjects who are contagious but not detected by that test.

**Author contributions** Conceptualization and methodology, GW, IT and RBW; Software and formal analysis, LL; Investigation and data curation, CWO, LH and RBW; resources, AG, GW, IT and RBW; Writing—original draft preparation, LL; Writing—review and editing, AG, CWO, GW, IT, LH and RBW; Supervision, GW. All authors have read and agreed to the final version of the manuscript.

**Funding** There was no funding source for this study.

## Compliance with ethical standards

**Conflict of interest** The authors declare no conflict of interest.

**Ethics approval** This study conformed to the principles outlined in the Declaration of Helsinki and was approved by ethics committee of the Innsbruck Medical University (ID of ethical vote: 1167/2020).

**Consent to participate** As the study was retrospective it was exempt from informed consent.

## References

1. World Health Organisation. Laboratory testing for coronavirus disease (COVID-19) in suspected human cases: WHO/COVID-19/laboratory/2020.5. <https://www.who.int/publications/i/item/10665-331501>
2. Bullard J, Dust K, Funk D, Strong JE, Alexander D, Garnett L, et al. Predicting infectious SARS-CoV-2 from diagnostic samples. *Clin Infect Dis*. 2020. <https://doi.org/10.1093/cid/ciaa638>.
3. Walsh KA, Jordan K, Clyne B, Rohde D, Drummond L, Byrne P, et al. SARS-CoV-2 detection, viral load and infectivity over the course of an infection. *J Infect*. 2020;81:357–71. <https://doi.org/10.1016/j.jinf.2020.06.067>.
4. Rao SN, Manissero D, Steele VR, Pareja J. A systematic review of the clinical utility of cycle threshold values in the context of COVID-19. *Infect Dis Ther*. 2020;9:573–86. <https://doi.org/10.1007/s40121-020-00324-3>.
5. Dzien A, Dzien-Bischinger C, Lechleitner M, Winner H, Weiss G. Will the COVID-19 pandemic slow down in the Northern

- hemisphere by the onset of summer? An epidemiological hypothesis. *Infection* 2020;48:627–9. <https://doi.org/10.1007/s15010-020-01460-1>.
6. Liotti FM, Menchinelli G, Lalle E, Palucci I, Marchetti S, Colavita F, et al. Performance of a novel diagnostic assay for rapid SARS-CoV-2 antigen detection in nasopharynx samples. *Clin Microbiol Infect.* 2020. <https://doi.org/10.1016/j.cmi.2020.09.030>.

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# Evaluation of the AESKU SARS-CoV-2 Antigen Rapid Test

## Purpose of the Study

The objective of this performance study is to establish the diagnostic sensitivity and diagnostic specificity of the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit (REF: 840001) and to provide data to demonstrate the product is effective for its intended use.

## Product Information

Manufacturer	AESKU.DIAGNOSTICS GmbH & Co. KG Mikroforum Ring 2 55234 Wendelsheim Germany Tel.: +49 6734 9622 0, <a href="mailto:info@aesku.com">info@aesku.com</a> , <a href="http://www.aesku.com">www.aesku.com</a>
Test Name	AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit (REF: 840001)
Detection Method	Immunochromatographic Test using a colored polymer-labeled novel coronavirus monoclonal antibody
Intended Use	Qualitative Detection of the N protein antigen from SARS-CoV-2 in human nasal swab specimen
Specimen	human nasal swab
Content of Testkit	AESKU.RAPID SARS-CoV-2 antigen test cassette Specimen processing tube Specimen sampling swab
Storage Condition	4-30°C
Lot number	P202010005 (exp date 09.04.2020), P202011003 (exp date 09.05.2020)

## Study Management

### *Sample Collection*

Biomex GmbH, Siemensstr. 38, 69123 Heidelberg, Germany

### *Investigation*

Hannes Deisel, Biomex GmbH, Siemensstr. 38, 69123 Heidelberg, Germany  
Sabine Szczesny, Biomex GmbH, Siemensstr. 38, 69123 Heidelberg, Germany  
Dr. Harald Depner, ESfEQA GmbH, Siemensstr. 38, 69123 Heidelberg, Germany

### *Study Coordinator and Author*

Dr. Dieter Groche, ESfEQA GmbH, Siemensstr. 38, 69123 Heidelberg, Germany  
Tel.: 06221 89466971, [groche@esfeqa.eu](mailto:groche@esfeqa.eu)

### *Testing Site*

Biomex GmbH, Siemensstr. 38, 69123 Heidelberg, Germany

### *Timelines*

The Study was performed in November and December 2020

## Study Design

### *Samples*

- 148 nasal swabs and 19 throat swabs from donors with known SARS-CoV-2 infection. Sex, age and symptoms of the donors as well as date of onset of symptoms were known. The date of infection was presumed from indications by the donor. Date of swab collections were documented. Samples were collected within 7 days after onset of symptoms.
- 164 nasal swabs and 50 throat swabs from healthy donors: Sex, age and date of sample collection were known (see annex "Comparison study SARS COV-2\_AESKU Rapid Test and PCR Test").

The nasal and throat swabs were collected between October 26th and December 12th and were stored at or below -20°C before they were analyzed.

### *Analytical investigation*

The study was performed on three separate days. Each sample was analyzed with the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test kit and immediately afterwards with the R-Biopharm (Darmstadt, Germany) real-time RT-PCR kit (see below). Identical sample preparations were used for both analytical methods.

Sample swabs were extracted in the AESKU.RAPID SARS-CoV-2 specimen processing tube as described in the IFU of the rapid test. Three drops of the specimen (approximately 145 µL) were added to the sample well of the test cassette. Results obtained with the rapid test device were visually read-out by two operators between 15 and 20 minutes after the sample had been applied onto the test cassette. Digital images were taken from used rapid test cassettes after visual read-out.

Total RNA was extracted from 300 µL of the remaining liquid using the R-Biopharm RIDA Xtract (REF:PGZ001) lot QL2000033 and lot QL200009, expiry date April 2022, and analyzed with the R-Biopharm RIDA Gene SARS-CoV-2 RUO real-time PCR kit (REF:PG6815 RUO) lot 24110N, expiry date March 2022, lot 26160N, expiry date April 2022, and RIDA Gene SARS-CoV-2 real time PCR kit (REF:PG6815) lot 24450N, expiry date November 2022. The instructions of the real-time RT-PCR kit manufacturer were followed with the exception that 300 µl instead of 400 µl of the solution was used for the extraction due to the limited volume in the specimen processing tube.

Real-time RT-PCR analysis was performed in duplicate analysis for all samples that were collected from infected donors and conducted using a CFX96 Touch Real-Time PCR Detection System from Bio-Rad Laboratories (Hercules, USA). The real-time RT-PCR results were obtained as Ct values. Samples with a Ct value above 36 (mean of the two replicates) were excluded from any statistical evaluation in this study.

## Results

In total 381 samples were tested in parallel with the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test and the R-Biopharm real-time RT-PCR assay and included in this study.

### Definitions

**True positive sample:** sample that was determined positive both using the AESKU.RAPID SARS-CoV-2 Antigen test and by RT-PCR.

**False positive sample:** sample that was determined positive using the AESKU.RAPID SARS-CoV-2 Antigen test, but negative by RT-PCR.

**True negative sample:** sample that was determined negative both using the AESKU.RAPID SARS-CoV-2 Antigen test and by RT-PCR.

**False negative sample:** sample that was determined negative using the AESKU.RAPID SARS-CoV-2 Antigen test but positive by RT-PCR.

**Specificity (%):**  $\# \text{ true negative samples} / (\# \text{ true negative samples} + \# \text{ false positive samples}) \times 100$

**Sensitivity (%):**  $\# \text{ true positive samples} / (\# \text{ true positive samples} + \# \text{ false negative samples}) \times 100$

Analytical Results for all samples with PCR result either negative or positive with a Ct value of less than 32:

		RT-PCR	
		positive	negative
AESKU.RAPID SARS-CoV-2 Antigen test	positive	105	4
	negative	4	218

Specificity of AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit: 98% (218/222), CI: 95-99%

Sensitivity of AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit (Ct < 32): 96% (105/109), CI: 91-99%

Analytical Results with correlation to Ct-values of the samples:

Ct value	Number of Samples	Number of true positive Rapid Test Samples	Number of false negative Rapid Test Samples	Sensitivity of AESKU.RAPID SARS-CoV-2 Antigen test (CI)
< 30	77	77	0	100 % (95-100)
< 32	105	101	4	96 % (91-99)
< 34	136	123	13	90 % (84-94)
< 36	157	133	24	85 % (78-90)

The correlation between the Ct-values of the analyzed samples and the sensitivity reveals a sensitivity of 100% for samples with a Ct-value of up to 30. Samples with a higher Ct value in the real-time RT-PCR and consequently less viral RNA copies as well as viral antigen in the samples result in lower

sensitivity values for the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test. This is in line with expectations regarding viral detection by antigen rapid testing compared to PCR analysis.

## Conclusion

The specificity and sensitivity of the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit was evaluated in this study with 381 samples collected as nasal or throat swabs. All samples were tested in parallel with the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit and a real-time RT-PCR assay. Samples with a Ct value below 32 were selected for the calculation of the sensitivity of the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit.

The specificity of the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit calculated from results of all samples was 98%, the sensitivity calculated from results of samples with a Ct-value less than 32 (105 samples) was 96% (95% CI: 91-99%). As expected, the sensitivity decreases by including samples with higher Ct value. Thus, by including all samples with a Ct value of or below 36 (157 samples) the sensitivity is calculated as 85% (95% CI: 78-90%).

In conclusion, the results from this study confirm that the AESKU RAPID SARS-CoV-2 Antigen Rapid Test Kit can be used for the qualitative detection of antigen from SARS-CoV-2 in human nasal swab and throat swab specimens.

## Approval

This study complements version 002 of the 'Evaluation of the AESKU SARS-CoV-2 Antigen Rapid Test' by including a larger number of samples.

Version 003 of the Evaluation Study Aesku Diagnostics SARS-CoV-2 Antigen Test was created and approved in December 2020 by the following persons:

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# Head-to-head comparison of SARS-CoV-2 antigen-detecting rapid test with self-collected anterior nasal swab versus professional-collected nasopharyngeal swab

Andreas K. Lindner<sup>1\*</sup>, Olga Nikolai<sup>1\*</sup>, Franka Kausch<sup>1</sup>, Mia Wintel<sup>1</sup>, Franziska Hommes<sup>1</sup>, Maximilian Gertler<sup>1</sup>, Lisa J. Krüger<sup>2</sup>, Mary Gaeddert<sup>2</sup>, Frank Tobian<sup>2</sup>, Federica Lainati<sup>2</sup>, Lisa Köppel<sup>2</sup>, Joachim Seybold<sup>3</sup>, Victor M. Corman<sup>4,5</sup>, Christian Drosten<sup>4,5</sup>, Jörg Hofmann<sup>6</sup>, Jilian A. Sacks<sup>7</sup>, Frank P. Mockenhaupt<sup>1,§</sup> and Claudia M. Denking<sup>2,§</sup>

<sup>\*,§</sup>Authors contributed equally

<sup>1</sup> Charité – Universitätsmedizin Berlin, corporate member of Freie Universität Berlin, Humboldt-Universität zu Berlin, and Berlin Institute of Health; Institute of Tropical Medicine and International Health, Berlin, Germany.

<sup>2</sup> Division of Clinical Tropical Medicine, Center of Infectious Diseases, Heidelberg University Hospital, Germany.

<sup>3</sup> Charité – Universitätsmedizin Berlin, corporate member of Freie Universität Berlin, Humboldt-Universität zu Berlin, and Berlin Institute of Health; Medical Directorate, Berlin, Germany.

<sup>4</sup> Charité – Universitätsmedizin Berlin, corporate member of Freie Universität Berlin, Humboldt-Universität zu Berlin, and Berlin Institute of Health; Institute of Virology, Berlin, Germany.

<sup>5</sup> German Centre for Infection Research (DZIF), partner site Charité, 10117 Berlin, Germany.

<sup>6</sup> Labor Berlin - Charité Vivantes GmbH, Berlin, Germany

<sup>7</sup> Foundation for Innovative New Diagnostics, Geneva, Switzerland

**Correspondence:** Claudia M. Denking, Division of Clinical Tropical Medicine, Heidelberg University Hospital, Im Neuenheimer Feld 672, 69120 Heidelberg, Germany.

E-mail: [claudia.denking@uni-heidelberg.de](mailto:claudia.denking@uni-heidelberg.de)

## Abstract

**Background:** Two antigen-detecting rapid diagnostic tests (Ag-RDTs) are now approved through the WHO Emergency Use Listing procedure and can be performed at the point-of-care. However, both tests use nasopharyngeal (NP) swab samples. NP swab samples must be collected by trained healthcare personnel with protective equipment and are frequently perceived as uncomfortable by patients.

**Methods:** This was a manufacturer-independent, prospective diagnostic accuracy study with comparison of a supervised, self-collected anterior nose (AN) swab sample with a professional-collected NP swab sample, using a WHO-listed SARS-CoV-2 Ag-RDT, STANDARD Q COVID-19 Ag Test (SD Biosensor), which is also being distributed by Roche. The reference standard was RT-PCR from an oro-/nasopharyngeal swab sample. Percent positive and negative agreement as well as sensitivity and specificity were calculated.

**Results:** Among the 289 participants, 39 (13.5%) tested positive for SARS-CoV-2 by RT-PCR. The positive percent agreement of the two different sampling techniques for the Ag-RDT was 90.6% (CI 75.8-96.8). The negative percent agreement was 99.2% (CI 97.2-99.8). The Ag-RDT with AN sampling showed a sensitivity of 74.4% (29/39 PCR positives detected; CI 58.9-85.4) and specificity of 99.2% (CI 97.1-99.8) compared to RT-PCR. The sensitivity with NP sampling was 79.5% (31/39 PCR positives detected; CI 64.5-89.2) and specificity was 99.6% (CI 97.8-100). In patients with high viral load ( $>7.0 \log_{10}$  RNA SARS-CoV2/swab), the sensitivity of the Ag-RDT with AN sampling was 96% and 100% with NP sampling.

**Conclusion:** Supervised self-sampling from the anterior nose is a reliable alternative to professional nasopharyngeal sampling using a WHO-listed SARS-CoV-2 Ag-RDT. Considering the ease-of-use of Ag-RDTs, self-sampling and potentially patient self-testing at home may be a future use case.

## *To the Editor:*

A number of antigen-detecting rapid diagnostic tests (Ag-RDTs) for SARS-CoV-2 are now commercially available and can result in rapid decisions on patient care, isolation and contact tracing at the point-of-care [1]. Two Ag-RDTs using nasopharyngeal (NP) swab samples meet WHO targets and are now approved through the WHO Emergency Use Listing procedure [2-4].

NP swab samples are frequently perceived as uncomfortable by patients and must be collected by trained healthcare personnel with protective equipment. Evidence supports the use of alternative sampling methods for RT-PCR, including anterior nasal (AN) swabs collected by patients and some tests have received regulatory approval with AN samples [5, 6]. Considering the ease-of-use of Ag-RDTs, a reliable simple sampling method would not only allow self-sampling, but may also pave the way for self-testing.

The primary objective of this prospective diagnostic accuracy study was a head-to-head comparison (positive and negative percent agreement) of a supervised, self-collected AN swab sample with a health care worker (professional)-collected NP swab sample, using a WHO-listed SARS-CoV-2 Ag-RDT against the reference standard RT-PCR collected from a NP/oropharyngeal (OP) swab. The secondary objective was to assess sensitivity and specificity for different sampling techniques with Ag-RDT. The study was continued until at least 30 positive NP swab samples according to Ag-RDT were obtained. This manufacturer-independent study was conducted in partnership with the Foundation of Innovative New Diagnostics (FIND), the WHO collaborating centre for COVID-19 diagnostics.

The study protocol was approved by the ethical review committee at Heidelberg University Hospital for the study site in Berlin (registration number S-180/2020). The study took place at the ambulatory SARS-CoV-2 testing facility of Charité University Hospital (Charité-Universitätsmedizin Berlin, Germany) from 23 September to 14 October 2020. The study enrolled adults at high risk for SARS-CoV-2 infection according to clinical suspicion. Participants were excluded if either of the swabs for the Ag-RDT or the RT-PCR reference standard could not be collected.

Participants underwent first an instructed, self-collected bilateral AN swab for the Ag-RDT. Verbal instruction was given to insert the swab horizontally 2-3 cm into the nostril and rotate it for 15 seconds against the nasal walls on each side. Deviations from the instructed technique were recorded. Subsequently, a combined OP/NP swab (eSwab from Copan) as per institutional recommendations for RT-PCR, and a separate NP swab for the Ag-RDT were taken from different sides of the nose. The samples for the Ag-RDTs were collected using the swab provided by the manufacturer within the test kit.

The Ag-RDT evaluated in this study was the STANDARD Q COVID-19 Ag Test (SD Biosensor, Inc. Gyeonggi-do, Korea; henceforth called STANDARD Q) [7], which is also being distributed by Roche [8]. The test uses the lateral flow assay principle in a cassette-based format with a visual read-out after 15-30 minutes. The manufacturer's instructions for use were followed. The Ag-RDTs were performed directly after sampling (within 60 minutes) at point-of-care by study physicians. The Ag-RDT results were interpreted by two operators, each blinded to the result of the other. The second reader was also blinded to the second Ag-RDT results of individual patients. The visual read out of the Ag-RDT test band was categorized on a semi-quantitative scale as negative, weak positive, positive and strong positive.

The Roche Cobas SARS-CoV-2 assay (Pleasanton, CA United States) or the SARS-CoV-2 E-gene assay from TibMolbiol (Berlin, Germany) were performed for RT-PCR according to routine procedures at the central laboratory. Viral RNA concentrations were calculated using assay specific CT-values, based on external calibrations curves [9, 10]. Staff performing the Ag-RDTs were blinded to results of RT-PCR tests and vice versa.

Of 303 patients invited, 289 (95.4%) consented to participate. Two patients were excluded as both swabs for the Ag-RDT could not be obtained. The average age of participants was 34.7 years (Standard Deviation [SD] 11.0) with 42.9% female and 19.0% having comorbidities. On the day of testing, 97.6% of participants had one or more symptoms consistent with COVID-19. Duration of symptoms at the time of presentation on average was 4.4 days (SD 2.7). Among the 289 participants, 39 (13.5%) tested positive for SARS-CoV-2 by RT-PCR (Table 1).

No invalid tests were observed on either AN or NP samples. Two patients were detected by NP but not by AN sampling. No patient was detected by AN sampling only. The positive percent agreement was 90.6% (CI 75.8-96.8; including 2 false positive results with AN and 1 with NP). The negative percent agreement was 99.2% (CI 97.2-99.8). Inter-rater reliability was near perfect with kappa of 0.98 for AN and 1.0 for NP samples. The semi-quantitative read-out was more often higher for the NP samples (9 higher on NP, 4 higher on AN). Of the two patients detected by NP but not by AN sampling, one patient collected the AN swab only with gentle rotation, the second presented 10 days post symptom onset with a low viral load (Table 1).

The STANDARD Q Ag-RDT with AN sampling showed a sensitivity of 74.4% (29/39 PCR positives detected; CI 58.9-85.4) and specificity of 99.2% (CI 97.1-99.8) compared to RT-PCR. The sensitivity with NP sampling was 79.5% (31/39 PCR positives detected; CI 64.5-89.2) and specificity was 99.6% (CI 97.8-100). In patients with high viral load ( $>7.0 \log_{10}$  RNA SARS-CoV2/swab), the sensitivity of the Ag-RDT with AN sampling was 96% (24/25 PCR positives detected; CI 80.5-99.8) and 100% (25/25 PCR positives detected; CI 86.7-100) with NP sampling. In contrast, the Ag-RDT frequently did not detect patients with lower viral load or with symptoms  $>7$  days (Table 1). For most patients, the application of the

flexible swab (meant for NP swab collection) in the anterior nose appeared unpleasant due to a tickling sensation and led to frequent sneezing.

The strengths of the study are the rigorous methods, including standardized sampling, two independent blinded readers and an additional semi-quantitative assessment of Ag-RDT results. The cohort was representative, judging from the comparable sensitivity observed in the recent independent validation study of STANDARD Q (sensitivity 76.6%; CI 62.8-86.4) [4]. The study is limited as it was performed in a single centre. Also, the NP swab was usually rotated against the nasopharyngeal wall for less time than recommended by the manufacturer, which may have a negative impact on the sensitivity of the Ag-RDT with NP sampling, but also reflects the difficulty of collection of this sample type.

In conclusion, this study demonstrates that supervised self-sampling from the anterior nose is a reliable alternative to professional nasopharyngeal sampling with STANDARD Q. The data will contribute to WHO recommendations for use of this test. Considering the ease-of-use of Ag-RDTs, self-sampling and potentially patient self-testing at home may be a future use case. If such testing could be repeated frequently and immediately ahead of situations when transmissions are likely to occur, self-testing with Ag-RDTs may have a significant impact on the pandemic. Further implementation studies on optimized self-sampling techniques and swabs (e.g. less flexible sponge swab) and the correct performance/interpretation of the test by patients themselves, are urgently needed to drive self-testing to scale.

**TABLE 1** Antigen-detecting RDT results with a supervised self-collected anterior nasal (AN) swab and with a professional-collected nasopharyngeal (NP) swab in RT-PCR positive patients from combined oro-/nasopharyngeal swab. CT-values and viral load (in descending order) of the paired RT-PCR samples are shown, as well as the duration of symptoms per patient. The positive percent agreement between AN and NP samples on Ag-RDT, as well as the respective sensitivities compared to RT-PCR are shown.

No.	AN swab self-collected SD Q Ag-RDT	NP swab prof.-collected SD Q Ag-RDT	OP/NP swab RT-PCR		Symptom duration (days)
			CT value	Viral load <sup>3</sup>	
1	pos (+++)	pos (+++)	17.33 <sup>1</sup>	9.92	2
2	pos (++)	pos (+++)	17.86 <sup>1</sup>	9.76	1
3	pos (+++)	pos (+++)	18.01 <sup>1</sup>	9.71	1
4	pos (++)	pos (+++)	18.31 <sup>1</sup>	9.62	3
5	pos (+++)	pos (+++)	18.40 <sup>1</sup>	9.60	3
6	pos (+++)	pos (+++)	18.76 <sup>1</sup>	9.49	4
7	pos (+++)	pos (+++)	18.77 <sup>1</sup>	9.49	5
8	pos (+++)	pos (+++)	18.78 <sup>1</sup>	9.49	5
9	pos (+++)	pos (+++)	19.05 <sup>1</sup>	9.41	3
10	pos. (+++)	pos. (+++)	19.40 <sup>1</sup>	9.30	2
11	neg.	pos (+++)	19.66 <sup>1</sup>	9.23	1
12	pos (+++)	pos (+++)	20.32 <sup>1</sup>	9.03	3
13	pos (++)	pos (+++)	20.44 <sup>1</sup>	9.00	2
14	pos (++)	pos (++)	20.54 <sup>1</sup>	8.96	5
15	pos (+++)	pos (+++)	17.81 <sup>2</sup>	8.86	4
16	pos (+++)	pos (+++)	21.09 <sup>1</sup>	8.80	4
17	pos (+++)	pos (+)	18.62 <sup>2</sup>	8.62	4
18	pos (+)	pos (++)	21.87 <sup>1</sup>	8.57	7
19	pos (++)	pos (+++)	22.05 <sup>1</sup>	8.52	2
20	pos (++)	pos (+++)	19.34 <sup>2</sup>	8.41	5
21	pos (+++)	pos (+++)	19.47 <sup>2</sup>	8.37	6
22	pos (+++)	pos (+++)	22.60 <sup>1</sup>	8.36	6
23	pos (+++)	pos (++)	23.66 <sup>1</sup>	8.04	6
24	pos (+)	pos (++)	26.42 <sup>1</sup>	7.23	5
25	pos (+++)	pos (+++)	26.77 <sup>1</sup>	7.12	5
26	neg.	neg.	24.25 <sup>2</sup>	6.96	10
27	pos (++)	pos (+++)	24.77 <sup>2</sup>	6.80	4
28	pos (+++)	pos (++)	25.29 <sup>2</sup>	6.64	2
29	pos (+)	pos (++)	29.33 <sup>1</sup>	6.36	5
30	neg.	neg.	29.56 <sup>1</sup>	6.30	3
31	neg.	neg.	29.95 <sup>1</sup>	6.18	3
32	pos (+)	pos (+)	30.25 <sup>1</sup>	6.09	4
33	neg.	neg.	27.81 <sup>2</sup>	5.90	8
34	pos (++)	pos (+)	31.20 <sup>1</sup>	5.81	8
35	neg.	pos (+)	31.61 <sup>1</sup>	5.69	10
36	neg.	neg.	32.58 <sup>1</sup>	5.40	10
37	neg.	neg.	32.86 <sup>1</sup>	5.32	2
38	neg.	neg.	34.62 <sup>1</sup>	4.80	7
39	neg.	neg.	35.53 <sup>1</sup>	4.53	14

Sensitivity  
29/39 (74.4%)

Positive percent agreement<sup>4</sup>  
90.6% (CI 75.8-96.8)

Sensitivity  
31/39 (79.5%)

<sup>1</sup> Roche Cobas SARS-CoV-2 assay (E-gene, T2 target)

<sup>2</sup> TibMolbiol assay, E-gene target.

<sup>3</sup> log<sub>10</sub> RNA SARS-CoV2/swab

<sup>4</sup> including 2 false positives on AN and 1 on NP

**Abbreviations:** No., patient number; SD Q, STANDARD Q COVID-19 Ag Test (SD Biosensor); Ag-RDT, antigen-detecting rapid diagnostic test; AN, anterior nasal; NP, nasopharyngeal; OP, oropharyngeal; CT, cycle threshold; RT-PCR, reverse transcription-polymerase chain reaction; neg., negative; pos (+), weak positive; pos. (++), positive; pos. (+++), strong positive.

**Acknowledgements:** Heike Rössig, Chiara Manon Rohardt, Claudia Hülso, Elisabeth Linzbach, Susen Burock, Katja von dem Busche, Stephanie Patberg, Melanie Bothmann, Zümrüt Tuncer, Stefanie Lunow, Beate Zimmer, Astrid Barrera Pesek, Sabrina Pein, Nicole Buchholz, Verena Haack, Oliver Deckwart.

**Author contributions:** AKL, LJK, FL and CMD designed the study and developed standard operating procedures. AKL and ON implemented the study design, enrolled patients, performed laboratory work and led the writing of the manuscript. FPM and JS coordinated and supervised the study site. FK, MW, FH enrolled patients. MGe coordinated the testing facility. MGa, LK and FT led the data analysis. VC, JH and CD were responsible for PCR testing and contributed to the interpretation of the data. JAS supported the study design setup and the interpretation of the data. CMD was the principle investigator of the study. All authors have reviewed the manuscript.

**Data availability:** All raw data and analysis code are available upon a request to the corresponding author.

**Conflict of interest:** None declared.

**Support statement:** The study was supported by FIND, Heidelberg University Hospital and Charité University Hospital internal funds, as well as a grant of the Ministry of Science, Research and the Arts of Baden-Württemberg, Germany. FIND provided input on the study design, and data analysis in collaboration with the rest of the study team.

## References

1. Foundation of Innovative New Diagnostics (FIND). SARS-CoV-2 Diagnostics Pipeline 2020 <https://www.finddx.org/covid-19/pipeline/>. Date last accessed: October 15 2020.
2. World Health Organisation. WHO interim guidance 11 September 2020. Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays 2020 [https://apps.who.int/iris/bitstream/handle/10665/334253/WHO-2019-nCoV-Antigen\\_Detection-2020.1-eng.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/334253/WHO-2019-nCoV-Antigen_Detection-2020.1-eng.pdf?sequence=1&isAllowed=y). Date last accessed: October 15 2020.
3. World Health Organisation. WHO news release 28 September 2020. Global partnership to make available 120 million affordable, quality COVID-19 rapid tests for low-and middle-income countries. <https://www.who.int/news-room/detail/28-09-2020-global-partnership-to-make-available-120-million-affordable-quality-covid-19-rapid-tests-for-low-and-middle-income-countries>. Date last accessed: October 15 2020.
4. Krueger LJ, Gaeddert M, Koeppl L, Bruemmer L, Gottschalk C, Miranda IB, et al. Evaluation of the accuracy, ease of use and limit of detection of novel, rapid, antigen-detecting point-of-care diagnostics for SARS-CoV-2. medRxiv. 2020:2020.10.01.20203836.
5. Tu YP, Jennings R, Hart B, Cangelosi GA, Wood RC, Wehber K, et al. Swabs Collected by Patients or Health Care Workers for SARS-CoV-2 Testing. N Engl J Med. 2020;383(5):494-6.
6. CDC. Interim Guidelines for Collecting H, and Testing Clinical Specimens for COVID-19. <https://www.cdc.gov/coronavirus/2019-ncov/lab/guidelines-clinical-specimens.html>. Date last accessed: October 15 2020.
7. SD Biosensor. COVID-19 Ag STANDARDTM Q COVID-19 Ag Test 2020. <http://sdbiosensor.com/xs/product/7672>. Date last accessed: October 15 2020.
8. Roche to launch SARS-CoV-2 Rapid Antigen Test in countries accepting CE mark, allowing fast triage decisions at point of care. <https://www.roche.com/media/releases/med-cor-2020-09-01b.htm>. Date last accessed: October 15 2020.

9. Matheeussen V, Corman VM, Donoso Mantke O, McCulloch E, Lammens C, Goossens H, et al. International external quality assessment for SARS-CoV-2 molecular detection and survey on clinical laboratory preparedness during the COVID-19 pandemic, April/May 2020. Euro Surveill. 2020;25(27).
10. Corman VM, Landt O, Kaiser M, Molenkamp R, Meijer A, Chu DK, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020;25(3).

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# HOW TO COLLECT YOUR ANTERIOR NASAL SWAB SAMPLE FOR COVID-19 TESTING

Follow the instructions included with your sample kit. Use **only** materials provided in your kit to collect and store your sample, unless the kit says to do otherwise. Use **only** an approved sample collection kit given to you by your healthcare provider or personnel at the testing center.

## Initial set-up

1. Open the sampling kit.

2. Apply hand sanitizer with at least 60% alcohol. Cover all surfaces of your hands and rub them together until they feel dry.

## Sample collection

3. Remove the swab from the container, being careful not to touch the soft end, which is the absorbent tip.

4. Insert the entire absorbent tip of the swab into your nostril, but do not insert the swab more than  $\frac{3}{4}$  of an inch (1.5 cm) into your nose.

5. Slowly rotate the swab in a circular path against the inside of your nostril at least 4 times for a total of 15 seconds. Be sure to collect any nasal drainage that may be present on the swab.

6. Gently remove the swab.

7. Using the same swab, repeat steps 4-6 in your other nostril.

## Preparation of sample for return

- 8.** Place the swab in the sterile tube and snap off the end of the swab at the break line. Place the cap on the tube.

- 9.** Re-apply hand sanitizer.

- 10.** Place the tube containing the swab in the biohazard bag provided and seal the bag.

## Returning the sample and clean-up

- 11.** Give the bag with the swab to testing personnel.




















- 12.** Throw away the remaining sample kit items.

- 13.** Re-apply hand sanitizer.

# Usability Study of the AESKU SARS-CoV-2 Antigen Rapid Test

## Purpose of the Study

The objective of this usability study is to evaluate the usability of the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit (REF: 840001) and to provide data to demonstrate that the product can be used safely by professional and laypersons.

## Product Information

Manufacturer	AESKU.DIAGNOSTICS GmbH & Co. KG Mikroforum Ring 2 55234 Wendelsheim Germany Tel.: +49 6734 9622 0, <a href="mailto:info@aesku.com">info@aesku.com</a> , <a href="http://www.aesku.com">www.aesku.com</a>
Test Name	AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit (REF: 840001)
Detection Method	Immunochromatographic Test using a colored polymer-labeled novel coronavirus monoclonal antibody
Intended Use	Qualitative Detection of the N protein antigen from SARS-CoV-2 in human nasal swab specimen
Specimen	human nasal swab
Content of Testkit	AESKU.RAPID SARS-CoV-2 antigen test cassette Specimen processing tube Specimen sampling swab
Storage Condition	4-30°C
Lot number	P202010005
Expiration Date	09.04.2022

## Study Management

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A questionnaire for the usability evaluation of the AESKU.RAPID SARS-CoV-2 Antigen test was developed by ESfEQA GmbH, Siemensstr. 38, 69123 Heidelberg, Germany

### Study Coordinator

ESfEQA GmbH, Siemensstr. 38. 69123 Heidelberg, Germany  
Tel.: 06221 89466971, [groche@esfeqa.eu](mailto:groche@esfeqa.eu)

### Timelines

The Study was performed in November 2020

## Study Design

### *Test persons*

Professionals experienced in IvD testing as well as laypersons were involved in this study.

### *Items provided to the test persons*

- AESKU.RAPID SARS-CoV-2 Antigen test, consisting of AESKU.RAPID SARS-CoV-2 antigen test cassette, specimen processing tube and specimen sampling swabs
- Gebrauchsanweisung AESKU.RAPID SARS-CoV-2 Rapid Test (in German language), Version 001 (Draft)
- Illustration of sample collection and test procedure (see annex)
- Questionnaire 'Fragebogen zum Usability-Test AESKU.Rapid SARS-CoV-2 Rapid Test' (in German language), Version 1, November 2020 (see annex)

### Results

Eleven individuals of the age between 16 and 84 years (5 female, 6 male) were selected for this usability study. The participants represent a wide range of age groups. Among them were 3 persons trained for performing in-vitro-diagnostic test. All participants were German-native speakers or at least fluent in German language.

The questionnaire addresses questions about the comprehension of the test procedure as well as the practical application of the test device. The individual questions (17 in total) could be answered in 5 different grades from very easy/perfect to very difficult/need of improvement.

All questions of the questionnaire were answered predominantly by the two highest grades (very easy/easy and perfect/almost perfect). The illustrations in the Instruction for Use ('Gebrauchsanweisung') and the description for the interpretation of test results were rated once as 'rather difficult'.

Several participants of this study (laypersons) missed the description of the sample collection in the chapter 'test execution' since they regard the sample collection as part of it. Instead, the sample collection is described in the IFU in the separate chapter 'sample material'. The separation of the description of the sample collection and testing procedure in the 1<sup>st</sup> version of the IFU may be reconsidered.

One participant of the study missed instructions in the IFU for the interpretation of a faint test line.

## Conclusion

The usability of the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit (REF: 840001) was evaluated in this study. The participants of this study, professionals in IvD testing as well as laypersons, have rated the comprehensibility as well as the practical test procedure predominantly as very easy to easy. Participants recommended to reconsider the separation of sample collection and test execution in two separate chapters in the IFU.

In conclusion, based on the predominantly positive feedback of the participants, the AESKU.RAPID SARS-CoV-2 Antigen Rapid Test Kit is evaluated as easy to use.

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